

Data Sheet:

Arial U™ Remote Push Button



Specifications	
Model	SKU: EAC-1000 (Cabletree)
Electrical	Max Voltage Rating: 30 VDC
Oxygen Safety	Meets UL 1069 Requirements
Latex Information	Latex Free
Mechanical	Button Life: 100,000 cycles Switch Type: Normally Open SP/ST
Size Measurements (Pendant)	Length: 3" Diameter: 1.25" (at widest point)

Listed to UL1069
File number: E496257
Type: HN



About STANLEY Healthcare

STANLEY Healthcare provides over 5,000 acute care hospitals and 12,000 long-term care organizations with enterprise solutions that create a safe, secure and efficient healthcare experience across life's stages. The STANLEY Healthcare solution set enables customers to achieve organizational excellence and superior care in critical areas: Patient/Resident Safety, Security & Protection, Environmental Monitoring, Clinical Operations & Workflow and Supply Chain & Asset Management. These solutions are complemented by STANLEY Healthcare's By Your Side™ Lifetime Customer Care commitment to ensure that every customer achieves success and realizes the full value of their investment, through consulting, training, implementation and integration services. STANLEY Healthcare is proud to be part of Stanley Black & Decker, Inc. For more information, visit stanleyhealthcare.com. Follow STANLEY Healthcare on [Facebook](#), [Twitter](#), [LinkedIn](#) and [YouTube](#).