


SOLUTION OVERVIEW

# WanderGuard® BLUE Wander Management

Helping caregivers protect residents

**WanderGuard®**  
BLUE



## Solution Benefits

- **Contemporary design**

Components are sleek and unobtrusive to blend into the décor. The resident wearable is one of the smallest on the market to preserve resident dignity.

- **Lower total cost of ownership**

Feature-rich components are wirelessly managed to simplify both physical installation and day-to-day management

- **Robust, scalable solution**

Designed for stand-alone exit protection or extended functionality and integrations to enable centralized alarm management and reporting

- **UL294-compliant access control**

Enable staff or authorized visitors to pass through a controlled area using special codes or access cards

- **Location-based alerts & notifications**

Easily identify and manage resident location. Geofencing creates virtual boundaries that can be used to initiate alerts based on a resident's presence and movement within the community.

## Challenges

Eloperments are among the most serious safety concerns for senior living communities. The challenge lies in providing protection for residents at risk of wandering while preserving residents' dignity and mobility. Memory care requires discreet, versatile technology that meets the safety and security needs of both residents and staff.

## Solution

Whether protecting one door or an entire campus, WanderGuard BLUE provides peace of mind, removing the worry of a resident wandering into danger. With a legacy of success, WanderGuard BLUE is the most sophisticated—yet most simple—way to protect memory care residents against elopement risk.

## How It Works


WanderGuard BLUE combines small and discreet resident bracelets and door components with streamlined design and simplified management. When a resident with a bracelet approaches a monitored door, the system alerts and the door locks when paired with optional magnetic door locks.

Advanced capabilities transform WanderGuard BLUE into a centralized alarm management and reporting solution, with the ability to manage daily wander management user activity from your desktop. Automatically identify resident location (based on tracking of supervised bracelets) with resident location clearly displayed on facility maps using the Arial® central management system. Geofencing creates virtual boundaries that can be used to initiate alerts based on a resident's presence and movement within the community.



*“STANLEY Healthcare’s integrated WanderGuard BLUE and Ariel solution empowers our caregivers with better tools to care for memory care residents. Providing location of residents and geofencing makes this wander protection solution ideal.”*

**MELISSA MILLER**  
**ASSISTED LIVING ADMINISTRATOR**  
**RIVERSIDE LODGE**



## About STANLEY Healthcare

Over 15,000 hospitals and senior living communities rely on STANLEY Healthcare solutions to empower caregivers to deliver better care. STANLEY Healthcare is a part of Stanley Black & Decker and a proud supporter of the Alzheimer’s Association®. Learn more at [stanleyhealthcare.com](https://stanleyhealthcare.com).